Lab 9. (case) 3D deconvolution of the microscopy fluorescence image

Materials

- Five '* mat' files
 - simSphere1.mat a non-blurred microsphere
 - PSFG&Lsmall.mat a calculated PSF data
 - simblurNoNoise.mat a calculated blurred microsphere without noise
 - Realdat.mat a measured microsphere data
 - PSFreal.mat a measured PSF data

Tasks

- Generate noise disturbed data using calculated blurred microsphere data;
 - Be careful, the original calculated blurred microsphere data assumes that the fluorescence intensity of each point is 1 and this fluorescence intensity should be adjusted to ensure the final value of image pixels are not exceed 1.
 - It is recommended to add background noise to simulate the real detector system.
- Implement the inverse filtering and Wiener filtering method;
 - MATLAB function 'deconvwnr' is not allowed.
- 3. Implement the R-L iteration algorithm. (Optional)

Tasks(cont.)

- 4. Compare the inverse filtering, Wiener filtering and R-L iteration on simulated data;
 - the R-L iteration method can applied based on either of TASK 3 and MATLAB function 'deconvlucy';
 - try to change the noise model and SNR to compare and discuss the different results.
- 5. Apply the inverse filtering, Wiener filtering and R-L iteration on real data and compare the results.

Methods and Tips

Inverse filtering

$$\overline{F}(u,v) = F(u,v) + \frac{N(u,v)}{H(u,v)}$$

$$\bar{f}(x, y, z) = \mathcal{F}^{-1} \left\{ \frac{G(\omega_x, \omega_y, \omega_z)}{H(\omega_x, \omega_y, \omega_z)} \right\}$$

Weiner filtering

$$\bar{F}(u,v) = \left[\frac{1}{H(u,v)} \frac{|H(u,v)|^2}{|H(u,v)|^2 + \frac{S_{\eta}(u,v)}{S_f(u,v)}} \right] G(u,v) = \left[\frac{1}{H(u,v)} \frac{|H(u,v)|^2}{|H(u,v)|^2 + K} \right] G(u,v)$$

The inverse filter in Weiner filtering

$$\widehat{H}(\omega_{x}, \omega_{y}, \omega_{z}) = \frac{H^{*}(\omega_{x}, \omega_{y}, \omega_{z})}{\left(\left|H(\omega_{x}, \omega_{y}, \omega_{z})\right|^{2} + \frac{S_{\eta}(\omega_{x}, \omega_{y}, \omega_{z})}{S_{f}(\omega_{x}, \omega_{y}, \omega_{z})}\right)}$$
NSR

Tips

•
$$S_f(\omega_x, \omega_y, \omega_z)$$

Power of the non-blurred signal

$$\frac{S_{\eta}(\omega_{x}, \omega_{y}, \omega_{z})}{S_{f}(\omega_{x}, \omega_{y}, \omega_{z})}$$

•
$$S_{\eta}(\omega_x, \omega_y, \omega_z)$$

Power of the noise

- When NSR = 0
 - Inverse filtering = Weiner filtering

Tips

- For simulated / calculated data
 - Be careful, the original calculated blurred microsphere data assumes that the fluorescence intensity of each point is 1 and this fluorescence intensity should be adjusted to ensure the final value of image pixels does not exceed 1.
 - It is recommended to add background noise to simulate the real detector system.
 - blurred * factor + background noise < 1
 - Power of (non-blurred signal*factor)
- For real data
 - Power of non-blurred signal?

Tips

$$\widehat{H}(\omega_{x}, \omega_{y}, \omega_{z}) = \frac{H^{*}(\omega_{x}, \omega_{y}, \omega_{z})}{\left(\left|H(\omega_{x}, \omega_{y}, \omega_{z})\right|^{2} + \frac{S_{\eta}(\omega_{x}, \omega_{y}, \omega_{z})}{S_{f}(\omega_{x}, \omega_{y}, \omega_{z})}\right)}$$

- Avoiding extremely small numbers in the denominator
 - Set these numbers to be numbers a little bigger
 - Useful Matlab functions: eps, realmin, ...
 - eps = 2.2204e-16
 - sqrt(eps) = 1.4901e-08
 - realmin = 2.2251e-308

Maximize the likelihood

Maximize the log likelihood of the Poisson process


$$\max J(f) = \sum_{x,y,z} g * log\lambda - \lambda$$

Apply the gradient-based iterative search algorithm(R-L iteration algorithm)


$$\bar{f}_{k+1} = \left[\frac{g}{\bar{f}_k \otimes h} \otimes (-h)\right] \bar{f}_k, \qquad (-h) = h(-x, -y, -z)$$

- imnoise
 - Add noise to image
- psf2otf
 - Convert point-spread function to optical transfer function
- fftn, ifftn
 - n-D fft, ifft
- flipdim
- deconvlucy


Examples



Examples



Examples



The end